



A comprehensive guide  
for selecting the right  
ZIPPIE® pediatric wheelchair.

## ZIPPIE® Pediatric WHEELCHAIR REFERENCE GUIDE

# ZIPPİE MATRIX

	ZIPPİE X'CAPE®	ZIPPİE® GS / GS SE	ZIPPİE ZONE®	ZIPPİE SPHYNX
--	----------------	--------------------	--------------	---------------

FRAME TYPE:	FOLDING	FOLDING or RIGID	RIGID	TILT-IN-SPACE Folding
-------------	---------	------------------	-------	-----------------------



Specifications are subject to change without notice.



KEY FEATURES:	<ul style="list-style-type: none"> <li>XLOCK® locking crossbrace for rigid performance and one-step push-to-lock and pull-to-fold features</li> <li>Z-finity® System for adjusting lower extremity length and angles</li> <li>VersaRail™ seat rail for shared placement of components</li> </ul>	<ul style="list-style-type: none"> <li>Built-in growth system</li> <li>Available in reverse configuration</li> <li>Proven design and the highest quality construction for trouble-free durability</li> </ul>	<ul style="list-style-type: none"> <li>Ultra lightweight</li> <li>Growth capabilities</li> <li>Easy adjustments</li> <li>Cool style</li> </ul>	<ul style="list-style-type: none"> <li>Adjustable static tilt settings of 10°, 20°, 30°</li> <li>Quick-adjust recline back with 8.5° to 100° range</li> <li>Ultra compact folded package for easy transport</li> </ul>
---------------	--	--	--	--

HCPCS CODE	E1238	Rigid: E1237 Folding: E1238	E1235, E1237	E1232, E1234, E1161
Transit Option	Yes <sup>1</sup>	Yes	Yes	Yes
Tilt Range	_____	_____	_____	10°, 20°, 30°
Chair Weight*	19 lbs. w/o footrests	Folding: 25 lbs. w/o footrests Rigid: 29 lbs. w/o footrest	14 lbs.	28 lbs. w/o footrests
Weight Capacity <sup>w/o seating</sup>	165 lbs.	165 lbs.	165 lbs.	250 lbs.
Overall Width*	15" to 26.5"	17" to 27"	15" to 25"	21" to 26"
Seat Width	8" to 16" (grows to 18")	10" to 18"	8" to 16"	10" to 18"
Seat Depth	8" to 18" (grows to 20")	10" to 20"	8" to 18"	10" to 19" (Includes 5" depth growth standard)
Back Angle Adjustment	Fore (-) 8° to Aft (+) 24°	Fore (-) 5° to Aft (+) 30°	Fore (-) 17° to Aft (+) 8°	85° to 100°
Front Seat-To-Floor Height	13" to 18"	14" to 20"	14" to 19"	16.5" to 21"
Rear Seat-To-Floor Height	11" to 18"	12.5" to 20"	13" to 19"	Tilt angle 10° = 16" Tilt angle 20° = 14" Tilt angle 30° = 11.5"
Footrest / Frame Angle	80° Swing-In/Out, 80° Fixed Front, Z-finity System: Adjustable ELR: 0° to 90°	Swing-In/Out: 60°, 70°, 80° HD Lift Off: 70°, 80°, 90° Swing-Away: 90° ELR, ALR: 0° to 90°	80°	Adjustable, removable: 0° to 120°
Back Type	Fixed Height, Height Adjustable, Dynamic	Fixed Height and Height Adjustable	Folding Angle Adjustable Double Lock-Down	Folding Angle Adjustable Double Lock-Down
Back Height	Fixed Height: 16", 20", 23" Ht Adjustable: 11" to 24", 12" to 24" Dynamic: 12" to 24"	Fixed Height: 19", 22", 25" Ht Adjustable: Rigid: 16" to 22", 19" to 25" Folding: 17" to 23", 20" to 26"	6" to 8" through 15" to 18"	Fixed Height: 24" or 26"
Center of Gravity Adjust.	3"	3"	2.5"	N/A
Rear Wheel Sizes	12", 16", 20", 22", 24"	12", 16", 18", 20", 22", 24"	18", 20", 22", 24"	12"
Caster Sizes	3", 4", 5", 6"	3", 4", 5", 6", 8"	3", 4", 5", 6"	6"
True Fit or One-Time Free Growth Kit	Yes	GS: Yes GS SE: No	No <sup>2</sup>	No
Warranty <sup>Refer to Owner's Manual for details</sup>	Frame: Lifetime Parts: 1 Year	Frame: Lifetime Parts: 1 Year	Frame: Lifetime Parts: 1 year	Frame: Lifetime Parts: 1 year

\* Actual measurements may vary according to the chair configuration and accessories.  
\*\* Without footrests and rear wheels

1. Conforms with ANSI/RESNA WC-19 standards. Please refer to Quickie® & Zippie® WC-19 Transit Securement Supplement (MK-100072) for additional information. Supplement is available for download at [www.sunrisemedical.com](http://www.sunrisemedical.com).

ZIPPIE® TS / TS SE	ZIPPIE IRIS®	ZIPPIE VOYAGE®
TILT-IN-SPACE Folding or Rigid	TILT-IN-SPACE Folding or Rigid	EARLY INTERVENTION Stroller



- |  |  |   |
|--|--|---|
| <ul style="list-style-type: none"> <li>• Tilt-in-space up to 45°</li> <li>• Lightweight welded aluminum frame construction ideal for portability</li> <li>• Designed to accept the widest possible variety of seating solutions</li> </ul> | <ul style="list-style-type: none"> <li>• Patented rotation-in-space technology ensures a smooth, easy-to-control tilting motion</li> <li>• Lightweight frame</li> <li>• Optional XLOCK® folding crossbrace for convenient transport</li> </ul> | <ul style="list-style-type: none"> <li>• Tilt-in-space range of 35°, optional 40° of recline with Advanced Seating and 90° of recline with Moderate Seating</li> <li>• Personalized seating with extensive built-in adjustments</li> <li>• Functional accessories for child and caregivers</li> </ul> |
|--|--|---|

Folding: E1234 Rigid: E1233	Rigid: E1233 Folding: E1234	E1232 without seating: E1234
Yes	Yes <sup>1</sup>	Yes
45°	40°: -5° to 35°, 0° to 40°, +5° to 45° 55°: -5° to 50°, 0° to 55°, +5° to 60°	35°: -5° to 30°
29 lbs. w/o footrests	Folding: 22.8 lbs. ** Rigid: 26.6 lbs. **	28 lbs.
165 lbs.	225 lbs. ◇	75 lbs.
17" to 26"	17.5" to 25.5"	24"
10" to 18"	Rigid: 10" to 18" Folding: 10" to 16" (grows to 18")	Adv.: 6" to 13" Mod.: 9.5" to 12"
13" to 20"	10" to 18"	Adv.: 5" to 13" Mod.: 8" to 13"
Fore (-) 5° to Aft (+) 30°	Fore (-) 5° to Aft (+) 30°	Adv.: 85° to 140° Mod.: 90° to 180°
15.5" to 20"	13" to 19" w/drop seat: 11" to 17.5"	_____
14" to 19"	13" to 19" w/drop seat: 11" to 17.5"	_____
Swing-In/Out: 60°, 70°, 80°, 90° HD Lift Off: 60°, 70°, 80° ELR, ALR: 0° to 90°	Swing-In/Out: 60°, 70°, 80°, 90° HD Lift Off: 70°, 80°, 90° ELR, ALR: 0° to 90° Contracture Z-finity System: Adjustable	0° to 90°
Fixed Height and Height Adjustable	Fixed Height, Height Adjustable, Reclining, Dynamic, MONO Backrest System®	Height Adjustable, Reclining
Fixed Height: 19", 22", 25" Ht Adjustable: Rigid: 16" to 22", 19" to 25" Folding: 17" to 23", 20" to 26"	Fixed Height: 19", 22", 25" Ht Adjustable: 16" to 22", 19" to 25" Reclining: 21", 24" MONO: 18.7", 21.7"	Adv.: Short: 7" to 13" Adv.: Tall: 10" to 17" Mod.: 16" to 24" (includes head support)
None	6"	_____
12", 16", 18", 20", 22", 24"	12", 16", 20", 22", 24"	12"
3", 4", 5", 6", 8"	4", 5", 6", 8"	7"
TS: Yes TS SE: No	IRIS: Yes IRIS SE: No	No
Frame: Lifetime Parts: 1 Year	Frame: Lifetime Parts: 1 Year	Frame: Lifetime Parts: 1 Year

2. Includes 3" of built-in frame depth growth and 2" of built-in seat sling depth growth.

◇ WARNING: Addition of the Z-Finity Footrest System to this device limits the maximum allowable user weight to 165 lbs. (75 kg)

ZIPPIE Q300 M Mini®
MID-WHEEL DRIVE POWER



- KEY FEATURES:**
- Our smallest power wheelchair base for superior maneuverability
  - Optional 12" power seat elevate, power tilt, and manual recline
  - Pediatric-specific seating with quick seat width and depth growth adjustments

HCPCS CODE	K0848 K0856
Transit Tested	Yes <sup>1</sup>
Drive Wheel	Mid
Max Speed	6.0 mph
User Weight Capacity	165 lbs.
Electronics	PG Drives VR2 or R-net
Suspension	6-wheel
Turning Radius	17"
Overall Base Width	20.5"
Overall Base Length	33" w/o footrests
Drive Wheels	12" or 14"
Casters	7"
Seat Width	10" to 16"
Seat Depth	10" to 19"
Seat-to-Floor Height	17" to 19"
Tilt	Power
Recline	Manual
Seat Elevate	12"
Batteries	Group 22

## ZIPPIE Frame Colors



### IMPORTANT NOTE:

The color options shown are for reference only. They are not intended as accurate representations.

Refer to the QUICKIE & ZIPPIE Color Brochure or a QUICKIE & ZIPPIE Color Ring for representations of ZIPPIE colors.

Print processes cannot accurately represent paint finishes. Actual colors may vary.

Colors and patterns may vary based on type of chair or frame material.

### Glow paint



## ZIPPIE Q300 M Mini Available Shroud Colors



## ZIPPIE Q300 M Mini Available Accent Colors

